

Nuclear Science and Security Consortium Virtual Scholar Showcase 2020

Organic Scintillator Characterization for Neutron Detection

June 3, 2020

Juan J. Manfredi University of California, Berkeley




June 2 - 3, 2020



Introduction









Sandia

National

Laboratories

Erik Brubaker

Bethany Goldblum

(Nuclear Engineering)

(California)


Optically Segmented Single-Volume Scatter Camera

- Scintillator characterization (UCB/LBNL)
- Simulation and imaging studies (SNL)

Focus Area: Radiation Detection and Instrumentation Crosscutting Area: Modeling and Simulation




- Organic scintillators can detect fast neutrons via np-elastic scattering
- Fast neutron detection central to NNSA Mission (SNM search/standoff detection, imaging, cargo screening, arms control, emergency response)
 o e.g., Single Volume Scatter Camera (Sandia, Brubaker)






- In this talk, I will discuss three different research efforts related to the characterization of organic scintillators
 - Proton light yield of several SVSC candidate materials
 - Novel organic glass from Sandia
 - Pulse height vs. pulse integral for proton light yield



PLY measurements at UCB/LBNL










- Several plastic scintillators from Eljen that are considered candidate materials for the SVSC were measured: EJ-200, EJ-204, EJ-208
- Results show similar light yields for all three (due to shared PVT base)
- Light yield was measured from 50 keV to 5 MeV, allowing for constraints on physics-based scintillation models







Proton light yields: EJ-23X




 Measurements of fast plastic materials also being considered for the SVSC (EJ-230, EJ-232, EJ-232Q) show proton light yield consistent with EJ-20x series (same PVT base)





- Does scintillator size affect the relative proton light yield?
- EJ-309 measurement in Enqvist, et al. NIM A 715 (2013) said yes
- Our work shows that there is no size dependence for EJ-232, EJ-232Q (which have short attenuation lengths)



- Novel melt-cast organic glass material has been developed at Sandia Livermore
- Relative electron light yield, proton light yield, and PSD have been characterized and compared to EJ-309 and EJ-276

J.S. Carlson and P.L. Feng, NIM A 832 (2016). T.A. Laplace, et al., IEEE Trans Nucl Sci (submitted).

Light Yield [relative to 477 keV e]

 10^{-2}

10

Novel organic glass developed at Sandia



Proton Energy [MeV]




Boron-loaded organic glass




• SNL glass can be loaded with boron to add thermal neutron detection via:

 $\label{eq:B} \begin{array}{ll} {}^{10}B+n = \\ \begin{cases} {}^{7}Li + \ {}^{4}He, & \mbox{Q= 2.792 MeV, 6\%} \\ {}^{7}Li + \ {}^{4}He + \gamma (477.6 \ \mbox{keV}), & \mbox{Q= 2.310 MeV, 94\%} \\ \end{cases}$

 Using 478 keV gamma + observation PSD, neutron scatter events and capture/gamma scatter events can be separated in target PSD (below)











- Light yields are sometimes reported using the pulse height (PH) of the digitized photodetector trace (or with a short integration length)
 - Strictly speaking, only the total pulse integral (PI) of a digitized trace is proportional to the number of scintillation photons
- Literature shows discrepancies in proton light yield for EJ-309
 - $\circ~$ Wide range of integration lengths used







- In organic scintillators pulse can be dependent on both particle and energy, implying that PH and PI are not proportional
- Comparison of proton light yields from PI and PH (or short integration length PI) show clear energy-dependent bias in relevant energy range
- Light yield measured with digitizer should use PI!








- Relative proton light yield for several organic scintillators
 - Consistent results across different materials with same PVT base
 - No size dependence found for EJ-232, EJ-232Q
- Novel organic glass from Sandia Livermore
 - Performs favorably (electron light yield, proton light yield, PSD) compared to commonly used PSD materials
 - Boron-loading allows for thermal neutron capture signal and outperforms commercially available EJ-254
- Pulse height vs. pulse integral
 - Proton light yields should be calculated using pulse integral, even for non-PSD materials



NSSC Experience + History









14

Acknowledgements



This material is based upon work supported by the Department of Energy National Nuclear Security Administration under Award Number DE-NA0003180.

Disclaimer: This presentation was prepared as an account of work sponsored by an agency of the United States Government. Neither the United States Government nor any agency thereof, nor any of their employees, makes any warranty, express or implied, or assumes any legal liability or responsibility for the accuracy, completeness, or usefulness of any information, apparatus, product, or process disclosed, or represents that its use would not infringe privately owned rights. Reference herein to any specific commercial product, process, or service by trade name, trademark, manufacturer, or otherwise does not necessarily constitute or imply its endorsement, recommendation, or favoring by the United States Government or any agency thereof. The views and opinions of authors expressed herein do not necessarily state or reflect those of the United States Government or any agency thereof.